STEEMA

DUAL TOT SLING SYSTEM

Instructions for use

ENGLISH

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FRANÇAIS

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DEUTSCH

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ESPAÑOL

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PORTUGUES

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ITALIANO

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TÜRKÇE

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INSTRUCTIONS FOR USE

DESCRIPTION

Steema is a kit indicated for the treatment of stress urinary incontinence, containing:

- 1 Steema Sling (Ref.: SL-100S), manufactured with biocompatible synthetic materials.
- 2 Deschamps-type needles (Ref.: DPN-ST), disposable, designed to be used together with the sling at the time of implantation.
- 1 Urethral Protector (Winged guide) (Ref.: GA) disposable, designed to be used optionally during the sling implantation.

All Steema components are provided sterile and ready-to-use.

The Steema Sling is a permanent implant consisting of a central polypropylene mesh covered by two protecting sheaths, with a positioning tab and two dilator arms with suture loops at their ends.

Once the surgical procedure has finished, the protecting sheaths, positioning tab, dilator arms and suture loops are removed. All of these are NOT implanted permanently.

INDICATION FOR USE

Steema has been designed for the treatment of female stress urinary incontinence due to urethral hypermobility and/or intrinsic sphincter deficiency.

CONTRAINDICATIONS

Steema should not be prescribed if there is any type of infection, especially genital or related to the urinary tract.

Steema should not be used in patients who:

- are pregnant,
- are sensitive or allergic to polypropylene products,
- have pre-existing pathologies or conditions posing an unacceptable surgical risk,
- have soft tissues pathologies in the site intended for implant placement, or
- have pathologies or other conditions that compromise wound healing.

WARNINGS

This product must be used only by properly trained surgeons with experience in the treatment of female stress urinary incontinence.

Careful patient screening, as well as a complete diagnostic study, is essential before surgery.

The surgical technique is safe, but in the event of severe deviations and/or major anatomical variations, there may occur perforations or injuries to blood vessels, organs or nerves, and surgical repair may be required.

When the needle passes, perforation or injury of blood vessels, organs or nerves may occur. Further repair surgery may be necessary.

Steema must be implanted without tension, that is, the central portion of the mesh must seat under the medium urethra without any tension. Implantation with tension brings about urinary retention and urethral erosion.

As with all foreign bodies, the polypropylene mesh could exacerbate an existing infection.

Steema components have been designed to be used only ONCE. Therefore, its components should NOT BE REUSED, REPROCESSED OR RESTERILIZED, as that may seriously harm the device performance and increase the risk of poor resterilization and cross contamination, which may result in patient injury, illness or death.

PATIENT INFORMATION

Surgeons are responsible for providing information to the patient or her representatives on possible complications associated to the sling implantation, prior to the surgery (see POSSIBLE COMPLICATIONS).

Promedon S.A. and their distributors delegate to the surgeon the responsibility to inform the patient of the advantages and possible risks related to the implantation and the use of the sling.

The patient must be warned that future pregnancies could invalidate the surgical effects of the sling implantation, and therefore, she might become incontinent again.

The patient should be warned that Steema is a permanent implant. Total or partial removal of the sling due to complications may imply additional surgeries.

It is recommended for the patient to avoid lifting weights and doing vigorous exercise involving effort (riding a bicycle, running, etc.), as well as maintaining sexual intercourse, during at least the first three or four weeks after surgery. The doctor should determine when it is appropriate for the patient to resume her normal activities.

The patient must immediately contact the surgeon in case of:

- Dysuria (pain or difficulty to urinate).
- Vaginal pain.
- Fever.
- Presence of serous, bloody or purulent secretion.

- Hemorrhages or other discomfort.
- Vaginal exteriorization of the mesh.

PRECAUTIONS

Steema must be used only by surgeons who are trained in the implantation technique and the postoperative management.

The sling must be handled with care, avoiding the use of pointed, serrated or sharp objects. Any damage, perforation or scratch may lead to complications. Fluff, fingerprints, powder, bacteria or other items contaminating the sling surface may cause infections or reactions to foreign bodies. Extra precautions must be taken to avoid contamination.

The surgical approach must be made carefully, avoiding major blood vessels and organs. Risks are minimized by paying attention to the local anatomy and passing the needle correctly.

The use of polypropylene in urogynecological procedures has been associated with cases of erosion, the treatment of which may require the partial surgical removal of the implant.

Like with all surgical procedures, there are certain known risk factors that may affect the results on the patient's pelvic floor and which may include, but are not limited to, altered vascularization (e.g. smoking status, estrogen status, previous radiation of the pelvic floor, etc.), or active infections in or near the surgical site. The patho-physiological conditions previously described must be considered at the time of determining whether the patient is a suitable candidate for sling implantation.

Hemorrhages may occur. The patient must be carefully monitored before leaving the hospital.

The operating room conditions must comply with the local hospital/administrative/government procedures and regulations.

After use, product and packaging must be disposed of pursuant to the local hospital/administrative/government procedures and regulations.

SUPPLY AND STORAGE

Steema is supplied sterile and pyrogen-free, in a pouch and/or blister.

Any damage to the sterile barriers renders the product non-sterile.

Storage must comply with the following conditions:

- TEMPERATURE: Room temperature.

DO NOT USE AFTER THE EXPIRATION DATE INDICATED IN PACKAGING.

MAGNETIC RESONANCE (MR) ENVIRONMENT

The Steema sling does not affect and is not affected by magnetic resonance (MR) environments.

POSSIBLE COMPLICATIONS

The possible complications associated with the use of the sling should be discussed with the patient prior to the surgery.

The use of this sling may result in complications related to the surgical procedure.

There may also exist complications associated with the patient's reaction or degree of intolerance to any foreign material implanted in her body.

Infections not responding to antibiotic treatment require the partial or total removal of the implant.

Some complications may require removal of the sling. Such procedure may involve multiple surgeries. In some cases, total removal of the sling may not be possible. Performance of multiple surgeries does not ensure the total repair of complications.

Some patients may experience groin or vaginal pain during the initial postoperative period. Treatment with ANALGESICS and ANTI-INFLAMMATORIES may be enough to relieve the pain.

Other complications reported with slings include, but are not limited to:

- Infection.
- Erosion of the urethral mucosa, vaginal mucosa, wall of bladder or other surrounding tissues.
- Pain or persistent pain (pelvic, vaginal, in the groin or thighs).
- Vaginal pain, discomfort or irritation.
- Dyspareunia.
- Purulent, serous or bloody discharge.
- Inflammation.
- Injury of blood vessels or nerves.
- Organ perforation or laceration.
- Bladder instability.
- Urinary obstruction.
- Urinary issues.
- Mesh exposure.
- Recurrence of incontinence.
- Hemorrhage.
- Neuromuscular issues.
- Vaginal healing.

Postoperative formation of a fibrous tissue capsule around the sling is a normal

physiological response to the implantation of a foreign body.

PROMEDON S.A. requires that all surgeons inform the Company or the Distributor of any complication observed with the use of Steema.

SURGICAL PROCEDURE

INSIDE-OUT APPROACH

The sling is implanted by a transobturator approach, following the currently accepted surgical technique, under regional, local or general anesthesia. The administration of prophylactic therapy with antibiotics should be considered, according to the procedure approved by the hospital.

The inside-out transobturator technique is summarized in the following steps:

1. Make a punctiform incision where a horizontal line passing at the level of the clitoris meets the genitofemoral folds.
2. Midline Colpotomy: Perform a sagittal incision, 1.5 cm long, starting about 1 cm from the lower edge of the urethral meatus.
3. Paraurethral Dissection: From the incision, perform a dissection at a 45° angle from the urethral axis, towards the obturator foramen in the paraurethral space. Take care not to injure the vaginal mucosa. Perform a minimal vaginal dissection in order to form a tunnel that will allow the passage of the transobturator needle. Repeat the process in both sides.
4. Thread the end of suture loop into the needle eyelet (see Figure 1).

Figure 1

5. Insert the urethral protector in the dissected area.
6. Enter the needle tip into the dissected area, following the urethral protector channel. Push the needle forward, slightly passing through and penetrating the obturator membrane.
7. Once in this position, remove the urethral protector and keep it sterile for later use in the same patient on the other side.
8. Rotate the instrument on the ischiopubic branch until the needle appears at the incision previously performed in the skin.
9. When the tip of the needle and the suture loop appear at the opening of the skin, unthread the loop, clamp it with forceps and remove the needle, rotating the handle in the opposite direction.
10. Pull the suture loop until it completely passes through the skin, and the sling with its protective cover appears.
11. Repeat the process in the other side of the patient (step 4 to 10). Make sure the positioning tab is located under the urethra and the sling is in horizontal position.
12. Cut the sling and its protective sheaths under the dilator arms. Locate the sling without tension under the mid-urethra by placing a Metzenbaum scissors between the sling and the urethra to make tension regulation easier and prevent the sling from bending. Pull the ends of the mesh until it makes contact with the scissors.
- When the sling is located without any tension under the urethra, remove the protective sheaths of the sling, without removing the Metzenbaum scissors.
13. Remove the Metzenbaum scissors.
14. Remove the positioning tab from the sling by cutting the thread located in the front. Make sure to take out the positioning tab and the thread from the vaginal channel.
15. Finally, cut the mesh excess and suture the incisions.

OUTSIDE-IN APPROACH

The sling is implanted by a transobturator approach, following the currently accepted surgical technique, under regional, local or general anesthesia. The administration of prophylactic therapy with antibiotics should be considered, according to the procedure approved by the hospital.

The transobturator technique is summarized in the following steps:

1. Make a punctiform incision where a horizontal line passing at the level of the clitoris meets the genitofemoral folds.
2. Midline Colpotomy: Perform a sagittal incision, 1.5 cm long, starting about 1 cm from the lower edge of the urethral meatus.
3. Paraurethral dissection: From the incision, perform a dissection at a 45° angle from the urethral axis, towards the obturator foramen in the paraurethral space. Take care not to injure the vaginal mucosa. Perform a minimal vaginal dissection in order to form a tunnel that will allow the passage

of the transobturator needle. Repeat the process in both sides.

4. Sling placement: The needle enters into the incision site towards the obturator foramen, bordering the ischiopubic branch until it reaches the vaginal incision. The surgeon should guide the exteriorization process by touching with his index finger to avoid urethral and vaginal mucosa injuries. Thread the suture loop through the needle eyelet (see Figure 2), and pull back the needle to pass the sling through the tunnel previously opened. Then, release the sling loop from the needle tip.

Figure 2

5. Repeat Step 4 in the other obturator hole.

6. Sling without tension: Place the Metzenbaum scissors between the urethra and the sling. Then, pull the suture loops, making sure the positioning tab is located under the urethra, preventing the sling from folding.

7. Cut the ends of the sling/sheaths under the dilator arms and remove the protective sheaths from both sides.

8. Remove the Metzembraum scissors.

9. Remove the positioning tab from the sling by cutting the thread located in the front. Make sure to take out the positioning tab and the thread from the vaginal channel.

10. Finally, cut the mesh excess and suture the incisions.

Postoperative care and therapy are at the surgeon's discretion.

In case a removal of implant is required, please note:

Polypropylene mesh integrate with patient's tissue, so complete removal may be difficult.

In case a mesh removal is necessary due to pain, we recommend trying to cut all the tension areas identified by the surgeon.

In most cases, the risk of organ injury caused by mesh removal may be higher than the benefits resulting from this removal, so each case should be assessed and decided at the surgeon's discretion.

SYMBOLS USED IN LABELS

	CATALOGUE NUMBER
	BATCH NUMBER
	SERIAL NUMBER
	CAUTION
	EXPIRY DATE
	DO NOT REUSE
	DATE OF MANUFACTURE
	STERILE. STERILIZATION METHOD: ETHYLENE OXIDE
	MANUFACTURER
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	DO NOT RESTERILIZE
	CONSULT INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGING IS DAMAGED
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY

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