Calistar

SINGLE INCISION POP REPAIR SYSTEM

Instructions for Use
ENGLISH

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Kullanim Talýmatlari
TÜRKÇE

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INSTRUCTIONS FOR USE

DESCRIPTION

Calistar A is a kit for the surgical treatment of the anterior and apical prolapse, that contains:

- 1 Calistar A implant (Ref.: CALISTAR-A), manufactured with biocompatible synthetic materials.
- 1 Retractable insertion guide (Ref.: DPN-MNL), disposable, designed to be utilized along with the implant when the placement is performed and provided with a tube in order to protect the tissues during the surgical procedure.
- 3 TAS (Tissue Anchoring System) anchors.

Calistar A consists of a permanent implant comprised of a central Type 1 single-filament polypropylene mesh with two polypropylene attachment arms and three TAS anchors with their corresponding sutures.

Calistar S is a kit for the surgical treatment of the anterior and apical prolapse that contains:

- 1 Calistar S implant (Ref.: CALISTAR-S), manufactured with biocompatible synthetic materials.
- 1 Retractable insertion guide (Ref.: DPN-MNL), disposable, designed for the placement of the anchors, and provided with a tube in order to protect the tissues during the surgical procedure.
- 1 Retractable insertion guide (Ref.: DPN-MNC) disposable, designed for the placement of the anterior attachment arms of the implant.
- 3 TAS (Tissue Anchoring System) anchors.
- 1 Knot pusher (Ref: KP), disposable, designed to be used during the surgical procedure if necessary, and provided with a protective tube.

Calistar S consists of a permanent implant comprised of a central Type 1 single-filament polypropylene mesh with two polypropylene attachment arms and three TAS anchors with their corresponding sutures. The central part of the implant is built with a low-weight, highly porous polypropylene mesh.

Calistar P is a kit for the surgical treatment of the posterior and apical prolapse, that contains:

- 1 Calistar P implant (Ref.: CALISTAR-P), manufactured with biocompatible synthetic materials.
- 1 Retractable insertion guide (Ref.: DPN-MNL), disposable, designed to be utilized along with the implant when the placement is performed and provided with a tube in order to protect the tissues during the surgical procedure.
- 3 TAS (Tissue Anchoring System) anchors.

Calistar P consists of a permanent implant comprised of a Type 1 polypropylene single-filament central mesh and three TAS anchors with their corresponding sutures.

All the components of Calistar A, Calistar S and Calistar P are supplied sterile and ready to use.

INDICATIONS

Calistar A and Calistar S

Calistar A and Calistar S are intended for the surgical treatment of the anterior and apical prolapse by means of tissue reinforcement and stabilization of the soft tissues of the female pelvic floor.

Calistar P

Calistar P is intended for the surgical treatment of the Posterior and Apical Prolapse by means of tissue reinforcement and stabilization of the soft tissues of the female pelvic floor.

CONTRAINDICATIONS

Calistar must not be used in patients

- undergoing anticoagulant therapy,
- with infectious processes, especially of the genital system or in the urinary tract,
- with vaginal, cervical, or uterine cancer,
- during pregnancy,
- with known sensitivity or allergy to polypropylene products,
- or with pre-existing conditions that pose an unacceptable surgical risk.

WARNINGS

Calistar must only be used by surgeons familiar with the procedures and techniques for implanting Calistar.

Careful patient selection and complete diagnostic evaluation are essential prior to surgery.

Maximum precautions must be taken to avoid contamination.

The surgical procedure must be carried out carefully, avoiding damage to large blood vessels, nerves, and organs. Risks can be minimized by closely observing the local anatomy and by inserting the retractable insertion guide correctly.

It has been established that the surgical technique is safe, but in the case of very severe deviations from this and/or very marked anatomic variations, when passing the retractable insertion guide, perforation or injury may occur to blood vessels, organs or nerves, and repair surgery may then be necessary.

As with all foreign bodies, the polypropylene mesh and attachment arms could exacerbate a pre-existing infection.

Calistar must be used with care in patients with:

- Diabetes.
- Coagulation problems.
- Obstruction of the upper urinary tract.
- Renal insufficiency.
- Autoimmune diseases affecting connective tissue.

Calistar components are supplied sterile. Do not use the product if the package is open or damaged.

Calistar components are designed for a SINGLE use only. Therefore, DO NOT REUSE or RESTERILIZE, since this could decrease the performance of the device and increase the risk of improper sterilization and cross-contamination.

The patient must be warned that future pregnancies could invalidate the surgical effects of the implant.

The patient should avoid heavy lifting and strenuous exercise that involves exertion (riding a bicycle, running, etc.) for the first three to four weeks after surgery, as well as sexual intercourse, for at least one month after the operation.

The patient should contact the surgeon immediately in the event of:

- Dysuria
- Vaginal pain
- Fever
- Serous, bloody or purulent secretions
- Hemorrhages or other problems

PRECAUTIONS

The implant must not be handled with pointed, serrated, or sharp objects since any damage, perforation, or tearing can cause subsequent complications.

Maximum precaution must be taken to avoid contamination.

Avoid excessive tension on the implant during insertion.

Maximum precautionary measures must be taken when joining the TAS to the insertion guide and to avoid exerting pressure in the wrong direction during TAS insertion. (See point "2" of the surgical procedure)

Operating room conditions must meet hospital, administrative or local government procedure.

After use, discard the product and packaging according to hospital, administrative or local government procedure.

SUPPLY AND STORAGE

Calistar is provided sterile and pyrogen-free in pouch or blister packaging. Any damage to the sterile barriers makes the device non-sterile.

Calistar should be stored under the following conditions:

- TEMPERATURE: Room temperature.

DO NOT USE AFTER THE EXPIRATION DATE INDICATED ON THE PACKAGING.

POSSIBLE COMPLICATIONS

Possible complications associated with the use of the implant should be discussed with the patient before surgery.

The use of a prosthesis may cause complications related to each patient's specific degree of intolerance to any foreign body implanted in the body. Some complications may require removal of the implant.

Some patients may experience vaginal pain during the initial post-operative period. Treatment with ANALGESICS and ANTI-INFLAMMATORY DRUGS may be sufficient to relieve pain.

Other reported complications with this or other similar implants include:

- Infection.
- Urethral or vaginal erosion.
- Mesh exposure.
- Adhesion formation.
- Vaginal pain, discomfort, irritation.
- Purulent, serous or bloody discharge.
- Inflammation.
- Injuries to blood vessels or nerves.
- Presence of vaginal fistula.
- Bladder instability.
- Urinary obstruction.
- Bowel problems.
- Recurrence of prolapse.

Post operative formation of fibrous tissue around the implant is a normal physiological response to the implantation of a foreign body.

In the event of vaginal exposure of the implant, generally due to infection, it may be necessary to partially remove it. This is achieved by cutting off the exposed part of the polypropylene mesh.

PROMEDON requires surgeons to report any complication associated with the use of Calistar to the Company or to the Distributor.

PATIENT INFORMATION

The surgeon is responsible for informing the patient or her representatives before surgery about the possible complications related to the implantation of Calistar.

The patient should be warned that future pregnancies could invalidate the surgical effects of implanting Calistar. The patient should avoid heavy lifting and

strenuous exercise involving exertion (riding a bicycle, running, etc.) and not have sexual relations until the doctor determines that it is appropriate to resume normal activities.

The patient should be advised that Calistar is a permanent implant, and any complication associated with the implant may or may not require additional surgery to correct the complication.

The patient should contact the surgeon immediately in the event of:

- Dysuria (pain or difficulty with urination).
- Vaginal pain.
- Fever.
- Serous, bloody or purulent secretions.
- Hemorrhages or other problems.
- Urinary obstruction.
- Bowel problems.

SURGICAL PROCEDURE

Prepare the patient for surgery in the usual manner and insert a Foley catheter into the urethra.

Calistar is implanted under local or regional anesthesia. Prophylactic antibiotic therapy should be administered according to procedures approved by the hospital.

The description of the technique is summarized in the following steps:

Calistar A

1. Urethrovesical Dissection:

Perform a dissection of the anterior vaginal wall. The incision of the anterior vaginal wall is made from the middle urethra to the cervix or the apex. The pubocervical fascia is carefully dissected. A blunt bilateral dissection is performed toward the ischial spine and coccygeal muscle, identifying the ischial spines, and then, the sacrospinous ligaments.

2. Insertion of the TAS (Tissue Anchoring System):

The TAS is inserted into the RIG (Retractable Insertion Guide).

Insert the TAS in the anterior wall of the sacrospinous ligament 2.5 cm medial to the ischial spine (Figure 1 shows the correct direction to apply pressure when inserting the TAS). The surgeon should use his index finger to touch and identify the ligament and to guide the RIG to its correct implant location. The TAS should be bilaterally placed, one in each sacrospinous ligament.


Figure 1

3. Implant placement

For the insertion of the implant in the suburethral zone, firstly, the RIG is connected to the front attachment arm and is introduced in the direction of the obturator internus muscle, guided by the surgeon's index finger. When the central mark of the implant is below the urethra, the retractable mechanism is activated to release the arm in its attachment place. This maneuver is repeated on the other side. A pair of Metzenbaum scissors can be introduced between the urethra and the implant in order to verify adequate suburethral tension and to avoid hypermobility overcorrection. Two attachment points are made with absorbable sutures on both sides of the mid urethra to prevent subsequent displacement of the mesh.

A cystoscopy is not required.

Then the polypropylene sutures of the TAS go through the reararms of the implant 1.5 cm from the end, taking care to keep a distance of at least 5 mm between both threads in order to make a double sliding knot.

Two additional attachment points must be made on the posterior body of the implant with non-absorbable sutures, which are fixed onto the remnants of the cardinal ligaments or the pericervical ring to prevent recurrence of cystocele.

Cystocele reduction is achieved by moving the implant towards the sacrospinous ligaments, using double knots previously made on the rear arms, to achieve full correction of the previous defect. The implant must remain below the bladder, free of tension and the rear end excess can be trimmed.

The vaginal incision is closed as usual. No cut is made in the vaginal wall unless necessary.

Calistar S

1. Bladder dissection:

Perform a dissection of the anterior vaginal wall. The incision of the anterior vaginal wall is made from the bladder neck towards the cervix or the apex. The pubocervical fascia is carefully dissected. A blunt bilateral dissection is performed toward the ischial spine and coccygeal muscle, identifying the ischial spines, and then, the sacrospinous ligaments.

2. Insertion of the TAS (Tissue Anchoring System):

The TAS is inserted into the corresponding RIG (Retractable Insertion Guide) (Ref.: DPN-MNL).

Insert the TAS in the anterior wall of the sacrospinous ligament 2.5 cm medial to the ischial spine (Figure 1 shows the correct direction to apply pressure when inserting the TAS). The surgeon should use his index finger to touch and identify the ligament and to guide the RIG to its correct implant location. The TAS should be bilaterally placed, one in each sacrospinous ligament.


Figure 1

3. Implant placement

For the insertion of the implant in the area of the bladder neck, first, the RIG for the anterior attachment arms (Ref.: DPN-MNC) is connected to the anterior attachment arm and is introduced towards the obturator internus muscle, the placement of the anterior attachment arms is made 1 cranial cm to the bladder neck at the balloon level. When the implant is in the correct position, the retractable mechanism is activated to release the arm in its attachment place. This maneuver is repeated on the other side. A pair of Metzenbaum scissors may be inserted between the bladder and the implant in order to ensure it is tension-free. Two attachment points are made with absorbable sutures on the central part of the implant to prevent subsequent displacement of the mesh.

A cystoscopy is not required.

Then the polypropylene sutures of the TAS go through the rear arms of the implant 1.5 cm from the end, taking care to keep a distance of at least 5 mm between both threads in order to make a double sliding knot.

Two additional attachment points must be made on the rear body of the implant with non-absorbable sutures, which are fixed in the remnants of the cardinal ligaments or the pericervical ring in order to favor apical suspension. Cystocele reduction is achieved by moving the implant towards the sacrospinous ligaments, using double nooses previously made on the rear arms, to achieve full correction of the previous defect. The implant must remain below the bladder, free of tension, and the rear end excess can be trimmed.

The vaginal incision is closed as usual.

Note: A plastic knot pusher and its protective tube are included in the kit in order to facilitate the knotting procedure of the TAS to the mesh, if the surgeon considers it convenient.

Calistar P

1. Pararectal Dissection

Perform a dissection of the posterior vaginal wall.

The incision in the posterior vaginal wall is made towards the cervix or the apex. Make a blunt bilateral dissection toward the ischial spine, then identify the coccygeal muscle and sacrospinous ligament on the right side. The same procedure is performed on the left side.

Sacrouterine ligaments are identified bilaterally and are held with Allis forceps. A 2-0 polypropylene suture point is applied to each sacrouterine ligament to reconstruct the posterior face of the pericervical ring.

2. Insertion of the TAS (Tissue Anchoring System):

The TAS is inserted into the RIG (Retractable Insertion Guide).

Insert the TAS in the anterior wall of the sacrospinous ligament 2.5 cm medial to the ischial spine (Figure 1 shows the correct direction to apply pressure when inserting the TAS). The surgeon should use his index finger to touch and identify the ligament and to guide the retractable insertion guide to its correct implant location. The TAS should be bilaterally placed, one in each sacrospinous ligament.


Figure 1

3. Implant placement

Then the polypropylene sutures of the TAS go through the rear arms of the implant 2 cm from the end, taking care to keep a distance of at least 5 mm between both threads in order to make a double sliding knot.

Two additional attachment points must be made on the posterior body of the

implant with non-absorbable sutures, which are attached to the sacrouterine ligaments and to the posterior surface of the pericervical ring. Reduction of the apical/posterior defect is achieved by moving the implant toward the sacrospinous ligaments, using the slipknots previously made on the arms of the implant, until the correction of the apical/posterior is achieved. The implant must remain below the bladder, free of tension, and the rear end excess may be trimmed. The front end excess is trimmed and then sutured bilaterally to the fascia of the levator ani muscle. The remnants of the rectovaginal fascia are used to cover the front part of the implant, which allows for greater protection from the exposure of the implant. The vaginal incision is closed as usual. Postoperative care and therapy are at the surgeon's discretion.

SYMBOLS USED IN LABELS

	CATALOG NUMBER
	LOT NUMBER
	SERIAL NUMBER
	CAUTION
	EXPIRY DATE
	DO NOT RE-USE
	MANUFACTURE DATE
	STERILE. STERILIZATION METHOD: ETHYLENE OXIDE
	MANUFACTURER
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
	CONSULT INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT RE-STERILIZE
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY

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