Unitape
t

Unitape
tplus

Unitape
VS

Instructions for use
ENGLISH

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Promedon

DESCRIPTION

Unitape t Plus, Unitape t and Unitape vs are kits for the treatment of stress urinary incontinence containing:

Unitape t Plus Kit

- 1 Unitape Sling (Ref: SL-100-L), manufactured with highly biocompatible synthetic materials.

- 2 Deschamps needles (Ref: DPN-HAS), single-use, designed to be used together with the sling for its implantation.

Unitape t Kit

- 1 Unitape Sling (Ref: SL-100-L), manufactured with highly biocompatible synthetic materials.

- 1 T needle (Ref: DPN-TAS), single-use, designed to be used together with the sling for its implantation.

Unitape vs Kit

- 1 Unitape Sling (Ref: SL-100-L), manufactured with highly biocompatible synthetic materials.

- 2 VS needle (Ref: DPN-SAS), single-use, designed to be used together with the sling for its implantation.

All Unitape t Plus, Unitape t and Unitape vs components are provided sterile and ready to use.

Unitape t Plus, Unitape t and Unitape vs slings are permanent implants that consist of a central polypropylene mesh that surrounds the urethra under the mid urethra, producing a coaptation of the former during strains.

INDICATIONS

The Unitape sling has been designed to treat female stress urinary incontinence caused by urethral hypermobility or intrinsic sphincter deficiency.

CONTRAINDICATIONS

Unitape cannot be prescribed if there is any kind of infection, specially genital or related to the urinary tract.

Unitape must not be used in patients

- undergoing anticoagulant therapy,
- with an existing urinary tract infection,
- during pregnancy,
- with known sensitivity or allergy to polypropylene products,
- or with pre-existing conditions that pose an unacceptable surgical risk,
- Unitape must be used cautiously in diabetic patients.

WARNINGS

The implant must not be handled with pointed, serrated, or sharp objects since any collapse, damage, perforation, or tearing can cause subsequent complications.

Fluff, fingerprints, talcum powder, or other contaminated surfaces may cause reactions to foreign bodies. Maximum precautions must be taken to avoid contamination.

Unitape must be implanted without tension, that is, the central mesh should rest on the urethra without any tension. An implantation with tension leads to urethral erosion and urinary retention problems.

Perforation or injury to blood vessels, organs or nerves could occur when passing the needle and another surgical repair may be needed.

As with all foreign bodies, the polypropylene mesh could exacerbate a pre-existing infection.

Unitape components are designed for SINGLE use only. Therefore, DO NOT REUSE or RESTERILIZE, since this could reduce the performance of the device and increase the risk of improper sterilization and cross-contamination.

PATIENT INFORMATION

The surgeon is responsible for informing the patient or her representatives before surgery about the possible complications related to the implantation of the sling.

Promedon S.A. and its distributors delegate to the surgeon the responsibility of informing the patient about the advantages and the possible risks related to the implantation and use of the sling.

The patient should be advised that the implantation of the sling might lose effectiveness with future pregnancies and that she could become incontinent again. The patient should avoid heavy lifting and strenuous exercise that involves exertion (riding a bicycle, running, etc.) for the first three to four weeks after surgery, as well as sexual intercourse, until at least one month after the operation.

The patient should contact the surgeon immediately in the event of:

- Dysuria (difficult or painful urination).
- Vaginal pain.
- Fever.
- Serous, bloody or purulent discharges.
- Hemorrhages or other difficulties.
- Vaginal exteriorization of the mesh.

PRECAUTIONS

Careful patient selection and complete diagnostic evaluation are essential prior to surgery.

The sling should be carefully manipulated avoiding the use of pointed, serrated or sharp objects.

Fluff, fingerprints, talcum powder, bacteria or other elements that contaminate the surface of the sling may cause infections or reactions to foreign bodies. Maximum precautions must be taken to avoid contamination. For that reason, it is recommended that the sling be kept in contact with antibiotics during surgery.

For that reason, it is specially recommended that the implant be placed in a sterile paste made of 10 gms of xylocaine jelly and two gentamicin ampoules (160 mg/ampoule). This mixture forms a gelatinous paste that adheres to the sling and when implanting it part of the mixture enters the body.

The surgical procedure must be carried out carefully, avoiding large blood vessels and organs. Risks are minimized if attention is paid to the local anatomy and if the needle is passed correctly.

Handling and Storage Precautions: Unitape is provided sterile and pyrogen-free. The packaging consists of a cardboard box that contains two pouches: one with the Sling and the other with the surgical instrument. IF ANY OF THE POUCHES ARE DAMAGED, DO NOT IMPLANT THE SLING.

Operating Room Conditions: In general, all operating rooms meet the requirements for the implantation of the sling. However, the following parameters should be especially examined:

- Asepsis of the Operating Room.
- Proper training of the operating room staff.

STORAGE

Unitape storage should meet the following conditions:

- TEMPERATURE: Room temperature

DO NOT USE AFTER THE EXPIRY DATE SPECIFIED IN THE PACKAGING.

POSSIBLE COMPLICATIONS

Possible complications associated with the use of sling should be discussed with the patient before surgery.

The use of the sling may result in complications associated with the medication and methods used during the surgical procedure, in addition to complications related to the patient's reaction or level of intolerance to any foreign body implanted in the body. Some complications may require removal of the sling. Infections that don't respond to antibiotic treatment require the partial or total removal of the prosthesis.

Some patients may experience suprapubic or vaginal pain during the initial post-operative period. Treatment with ANALGESICS and ANTI-INFLAMMATORY DRUGS may be sufficient to relieve pain.

Other complications reported with this or other slings include:

- Wound infection.
- Urethral or vaginal erosion.
- Vaginal pain.
- Purulent, serous or bloody discharge.
- Vaginal inflammation.
- Injuries to blood vessels or nerves.
- Presence of vaginal fistula.
- Bladder instability.
- Vaginal dehiscence.
- Urinary obstruction.

The postoperative formation of a fibrous tissue capsule around the sling is a normal physiological response to the implantation of a foreign body.

Promedon requires surgeons to report any complication associated with the use of Unitape to the Company or to the Distributor.

SURGICAL PROCEDURE

Prepare the patient for surgery in the usual manner and insert a Foley catheter into the urethra.

Unitape t Plus

The sling is implanted with a transobturator approach, under regional or local anesthesia, following the currently accepted surgery technique. The use of intravenous prophylactic antibiotics should be considered.

The description of the transobturator technique is summarized in the following steps:

1. Median Colpotomy

Make a sagittal incision, 1.5 cm long, starting approximately 1 cm from the outer urethral meatus.

2. Paraurethral Dissection

From the incision, release the vaginal wall with scissors. Make a minimal vaginal dissection in order to form a tunnel that allows the passage of the transobturator needle.

3. Sling Placement

Make a punctiform cutaneous incision at clitoris level, on the same horizontal

plane, in the genitofemoral folds on both sides.

Needles are introduced parallel to the ischiopubic ramus and exteriorized through the vaginal incision, simply by rotation of the fists, guided by the surgeons forefinger.

Thread the end of the sling through the tip of the needle and transfer the connector through the tunnel previously created with the needle.

Repeat this step in the other obturator orifice.

4. Tension-free fixation

Place a pair of Metzenbaum scissors between the mesh and the urethra to facilitate tension regulation and to prevent the mesh folding. Pull the connectors upwards until the mesh comes into contact with the urethra.

The correct tension can be determined by having the patient cough to verify continence.

Finally, cut the mesh excess and suture the incisions.

Postoperative care and therapy are at the surgeon's discretion.

Unitape t

The sling is implanted with a transobturator approach, under regional or local anesthesia, following the currently accepted surgery technique. The use of intravenous prophylactic antibiotics should be considered.

The description of the transobturator technique is summarized in the following steps:

1. MEDIAN COLPOTOMY

Make a sagittal incision, 1.5 cm long, starting approximately 1 cm from the outer urethral meatus.

2. PARAURETHRAL DISSECTION

From the incision, release the vaginal wall with scissors. Make a minimal vaginal dissection in order to form a tunnel that allows the passage of the transobturator needle.

3. SLING PLACEMENT

Make a punctiform cutaneous incision at clitoris level, on the same horizontal plane, in the genitofemoral folds on both sides.

Introduce the transobturator needle vertically through the cutaneous incision, passing through the internal transobturator membrane and the muscle, with the tip downward to avoid vascular injuries. Bring the needle to a horizontal position, behind the ischiopubic ramus, guide it towards the urethra, keeping a finger in the vaginal incision to protect the vaginal wall from the tip of the needle. Remove the needle through the vagina.

Thread the end of the sling through the tip of the needle and transfer the sling through the tunnel previously created with the needle.

Repeat this step in the other obturator orifice.

4. TENSION-FREE FIXATION

Place a pair of Metzenbaum scissors between the mesh and the urethra to facilitate tension regulation and to prevent the mesh folding. Pull the ends of the sling upwards until the mesh comes into contact with the urethra.

The correct tension can be determined by having the patient cough to verify continence.

Finally, cut the mesh excess and suture the incisions.

Postoperative care and therapy are at the surgeon's discretion.

Unitape vs

The sling can be implanted either with a vaginal or suprapubic approach, following any of the currently accepted and used surgical techniques, under regional or local anesthesia. The use of intravenous prophylactic antibiotics should be considered.

A. VAGINAL APPROACH STEPS (BOTTOM-UP)

1. MEDIAN COLPOTOMY

Make a sagittal incision, 1.5 cm long, starting approximately 1 cm below the urethral meatus.

2. PARAURETHRAL DISSECTION

Paraurethral dissection into the retropubic space, through sharp and blunt dissection.

3. ASSEMBLE NEEDLES WITH HANDLES

Introduce the hooked tip of the needle into the handle, leaving the conventional tip of the needle free.

4. Use your index finger to move the mid urethra and protect it from possible injury during the needle insertion through the vagina towards the suprapubic area.

5. Introduce the needle through the vaginal incision into the retropubic space, keeping the needle tip in constant contact with the posterior surface of the pubis and towards the patient's ipsilateral shoulder, until the free end of the needle appears in the suprapubic punctures. These are located bilaterally 2-3 cm from the midline.

6. The same procedure is performed on the other side of the patient.

7. Perform cystoscopy to ensure the bladder has not been perforated.

8. Remove the handles and connect them again to the conventional ends of the needles this time.

9. Thread the end of the sling at the tip of the needle and pull it up to transfer the sling to the suprapubic area.

10. Repeat the same maneuver on the other side of the patient.

11. TENSION-FREE FIXATION

Place Metzenbaum scissors between the sling and the urethra to ensure it is tension-free and prevent the sling from folding in on itself. Pull both arms of the sling up until it makes contact with the urethra. The correct tension can be determined by carrying out a Cough Test to verify continence.

12. Cut excess sling arms and suture the incisions.

B. SUPRAPUBIC APPROACH STEPS (TOP-BOTTOM)

1. Make a sagittal incision, 1.5 cm long, starting approximately 1 cm below the urethral meatus.

2. PARAURETHRAL DISSECTION

Paraurethral dissection into the retropubic space, through sharp and blunt dissection.

3. ASSEMBLE NEEDLES WITH HANDLES

Introduce the conventional end of the needle into the handle, leaving the hooked end of the needle free.

4. PUNCTURE IN SUPRAPUBIC AREA

Perform two punctures in the skin with a scalpel over the upper edge of the pubis, between 2 to 3 cm on each side of the midline. Separation between punctures: approximately 5 to 6 centimeters.

Introduce the needle in the suprapubic puncture towards the vaginal incision, always very close to the posterior wall of the pubis. The needle will make contact and will be guided by the surgeon's finger, previously introduced into the vaginal incision.

5. The same procedure is performed on the other side of the patient.

6. Perform cystoscopy to ensure the bladder has not been perforated.

7. Thread the end of the sling at the hooked needle tip and pull it up to transfer the sling to the suprapubic area.

8. Repeat the same maneuver on the other side of the patient.

9. TENSION-FREE FIXATION

Place Metzenbaum scissors between the sling and the urethra to ensure it is tension-free and prevent the sling from folding in on itself. Pull both arms of the sling up until it makes contact with the urethra. The correct tension can be determined by carrying out a Cough Test to verify continence.

12. Cut excess sling arms and suture the incisions.

Attention:

Control cystoscopy is recommended. In the case of vesical perforation during the suprapubic puncture, it is recommended to puncture the affected side again, further out from the first puncture. This is to create a new transfer path for the sling, leaving enough tissue between the sling and the punctured area of the bladder to enable the orifice to close and for the body to spontaneously protect itself better against any possible infection. In this case, leave a urethral probe for 10 to 15 days.

Postoperative care and therapy are at the surgeon's discretion.

SYMBOLS USED IN LABELS

PRODUCT CODE

LOT NUMBER

CAUTION

EXPIRY DATE

DO NOT REUSE

DATE OF MANUFACTURE


STERILE. STERILIZATION METHOD:
ETHYLENE OXIDE

DO NOT RESTERILIZE

MANUFACTURER

AUTHORIZED REPRESENTATIVE IN
THE EUROPEAN COMMUNITY

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