


# BIOMESH P1

Non resorbable parietal reinforcement implant




---

**Best** of the classics

---

# BIOMESH P1 : NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT

## > KEY POINTS



- Resistant
- Best of the classics

## > REFERENCES AND SIZES

	Shape	Reference	Description and size (cm)	
Lichtenstein		FBIO P1 6011	Pre-cut mesh with slit 6 x 11	
		FBIO P1 6135	Pre-cut mesh with slit 6 x 13	
		FBIO P1 7515	Pre-cut mesh with slit 7 x 15	
		FBIO P1 T510	Pre-cut mesh with slit 5 x 10	
TAPP/TEP		FBIO P1 0510	Mesh 5 x 10	
		FBIO P1 0611	Mesh 6 x 11	
		FBIO P1 0715	Mesh 7 x 15	
		VCB P1 6011X	Pre-cut mesh with slit 6 x 11 (10/box)	
		FBIO P1 1015	Mesh 10 x 15	
		FBIO P1 1515	Mesh 15 x 15	
		VCB P1 1515X	Mesh 15 x 15 (10/box)	
		XBIO P1 1515	Mesh 15 x 15 (10/box)	
	Ventral (extra-peritoneal)		FBIO P1 1015	Mesh 10 x 15
			FBIO P1 1515	Mesh 15 x 15
FBIO P1 1530*			Mesh 15 x 30	
FBIO P1 2536			Mesh 25 x 36	
		FBIO P1 3030	Mesh 30 x 30	
		FBIO P1 4545*	Mesh 45 x 45	
		VBIO P1 3030	Mesh 30 x 30 (5/box)	
		VCB P1 1515X	Mesh 15 x 15 (10/box)	
	XBIO P1 1515	Mesh 15 x 15 (10/box)		

\* : Manufactured on request

## > MATERIAL AND WEIGHT

Polypropylene  
100 ± 15 g/m<sup>2</sup>

Biomesh P1 is a class IIb medical device manufactured by COUSIN BIOTECH s.a.s. The management system of COUSIN BIOTECH s.a.s has been certified as meeting the requirements of ISO 13485. Please read instructions for use carefully. Reference: FPVMP1GB01 - 24/02/17

Non contractual pictures and texts. Specifications likely to be modified without notice.

Cousin Biotech S.A.S capital : 340 656 € - 398 460 261 RCS Lille - N°TVA FR 34 398 460 261

CE  
0120

